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
1. Ronald E. Wrolstad, Ph.D., "The Possible Health Benefits of Anthocyanin Pigments and Polyphenolics," The Linus Pauling Institute Oregon State University, <http://lpi.oregonstate.edu/ss01/anthocyanin.html>.
2. J. P. Goiffon, M. Brun, M. J. Bourrier, J. Chromatography. 537,101-121, 1991.
3. Robert D. Ricker, "High-Efficiency and High-Speed Separation of Natural Anthocyanins," Agilent Technologies, April 25, 2002, 5988-6362EN.
4. W. Kalt, J. E. McDonald, R. D. Ricker, and X. Lu, "Anthocyanin content and profile within and among blueberry species," Can. J. Plant Sci. 79: 617-623.
5. Anne E. Mack and William J. Long and "Rapid Tea Analysis on Poroshell 120 SB-C18," April 14, 2011, Agilent Technologies publication 5990-7842EN.
6. Anne E. Mack and William J. Long "Fast, Low Pressure analysis of food and beverage additives using a Superficially Porous Agilent Poroshell 120 EC-C18 Column," July 8, 2010, Agilent Technologies publication 5990-6082EN.
7. Agilent Technical Note, "Transfer of Methods between Poroshell 120 EC-C18 and ZORBAX Eclipse Plus C18 Columns," February 9, 2011, Agilent Technologies publication 5990-6588EN.
8. Dorothy J. Phillips, Mark Capparella, Uwe D. Neue, Zoubair El Fallah, "A new small particle packing for faster analysis with high resolution," Journal of Pharmaceutical and Biomedical Analysis, 15 (1997) 1389-1395.
9. William Long and Anne E. Mack, "Fast Analysis of Environmental Phenols with Poroshell 120 EC-C18 Columns," August 31, 2010, Agilent Technologies publication 5990-6156.
10. Xianli Wu and Ronald L. Prior, Systematic Identification and Characterization of Anthocyanins by HPLC-ESI-MS/MS in Common Foods in the United States: Fruits and Berries, J. Agric. Food Chem. 2005, 53, 2589-2599.
11. Ronald W. Durst and Ronald E. Wrolstad, "Separation and Characterization of Anthocyanins by HPLC," Current Protocols in Food Analytical Chemistry
12. Xianli Wu, Gary R. Beecher, Joanne M. Holden, David B. Haytowitz, Susan E. Gebhardt, and Ronald L. Prior, "Concentrations of Anthocyanins in Common Foods in the United States and Estimation of Normal Consumption", J. Agric. Food Chem. 2006, 54, 4069-4075.
13. Victor Hong and Ronald E. Wrolstad, "Use of HPLC Separation/Photodiode Array Detection for Characterization of Anthocyanins," J. Agric. Food Chem. 1990, 38, 708-715.

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© Agilent Technologies, Inc., 2011